



DIY Tape Microscope Slides Activity

Crystal Comparisons: Is it all that sweet?

Materials:

- Clear Tape
- Salt
- Sugar table /raw
- Artificial Sweetener(s): Splenda, Sweet and Low, Equal, Stevia
- Felt Tip Marker

Activity Steps:

1. Tear off a piece of clear tape
2. Fold over both sides of the tape to make "handles"
3. Place sample(s) on the adhesive area of the slide
4. Label your sample(s)
5. View under the microscope or through your tablet using the Moticam X (on a microscope) and MotiConnect App

Make your observations:

1.

2.

3.

Conclusions:

1.

2.

3.

